



MEDITECH
INDIA PVT. LTD.

over
30
YEARS
SERVICE & EXPERIENCE



EL – 580

Technical Specification

Light head dia: 15 cm

Illumination intensity: 40,000 Lux $\pm$ 10%

Focus dia: 15-20 cm

Color temperature: 4500 K

CRI (Colour Rendering Index): 93

Brightness control: Touch Panel

No. of LEDs: 14

LEDs lifetime: 30000 hrs.

Power supply: 220V/50Hz AC

📍 #487/67, Glass Street, National Market,
Peeragarhi, New Delhi - 110087, INDIA

📞 +91 11 40392484

☎ +91 9812266676

✉ marketinguk@meditech-india.com

🌐 www.mediech-india.com

A trusted Indian manufacturer and global exporter of Operating Tables, Surgical & Examination Lighting, Suction Machine, Fogger and Hospital Furniture - serving healthcare since 1996