

C-Arm Compatible Hydraulic Mechanical Operated OT Table

MI 111



Option : Semi Automatic (electric height adjustment) version is available



FEATURES

- Hydraulic and manual drives, thus very smooth positioning with no oil leakage.
- Table positions i.e. Trendelenburg, Reverse Trendelenburg, Lateral Tilt, Back Section/Chair Position, Flex-Reflex are controlled by gear operated side end crank.
- Table top is c-arm & x-ray compatible over the entire length of the table.
- Detachable-Interchangeable head and leg sections provide further versatility.
- Up & Down positions is activated by foot pedal operated concealed hydraulic pump.
- Non-corrosive stainless steel 304G construction including table top frames, accessories rail, base cover and column covers.
- Four easy glide, antistatic castors, two of which are swivelling, ensure mobility.

Technical Constraints

Top Dimension	L 1980 mm x W 533 mm
Height Adjustment	750 mm - 1050 mm
Trendelenburg/ Reverse	30° / 25°
Lateral Tilt	20° / 20°
Kidney Elevator	150 mm
Back Rest (Up/Down)	80° / 25°
Leg Rest (Up/Down)	15° / 90°
Head Rest (Up/Down)	20° / 60°
Patient Weight Capacity	150 kgs(330lbs)

Standard Accessories

- Anaesthesia screen – 1 pc
- Lateral support – 1 pair
- Shoulder support – 1 pair
- Arm rest – 1 pair
- Knee crutches – 1 pair
- Flat clamp – 7 pc
- Radial clamp – 1 pair
- Mattress 50 mm – 1 set