

Hydraulic-Mechanical Operated OT Table

MI 2003



Option : Semi Automatic (electric height adjustment) version is available



FEATURES

- Hydraulic and manual drives, thus very smooth positioning with no oil leakage.
- Table positions i.e. Trendelenburg, Reverse Trendelenburg, Lateral Tilt, Back Section/Chair Position, Flex-Reflex are controlled by gear operated side end crank.
- Detachable head and leg sections provide further versatility.
- Up & Down positions is activated by foot pedal operated concealed hydraulic pump.
- Base cover, column telescopic covers and accessories rails are made of Stainless steel 304G.
- Four easy glide, antistatic castors, two of which are swivelling, ensure mobility.

Technical Constraints

Top Dimension	L 1830 mm x W 497 mm
Height Adjustment	750 mm - 1000 mm
Trendelenburg/ Reverse	30° / 25°
Lateral Tilt	20° / 20°
Kidney Elevator	150 mm
Back Rest (Up/Down)	80° / 25°
Leg Rest (Up/Down)	15° / 90°
Head Rest (Up/Down)	20° / 60°
Patient Weight Capacity	120 kgs(264lbs)

Standard Accessories

- Anaesthesia screen – 1 pc
- Lateral support – 1 pair
- Shoulder support – 1 pair
- Arm rest – 1 pair
- Knee crutches – 1 pair
- Flat clamp – 7 pc
- Radial clamp – 1 pair
- Mattress 50 mm – 1 set