

Orthopaedic Attachment

MI 112



FEATURES

- The orthopaedic attachment can be used with any operating table.
- In addition, a pelvic support plate is supported on top of the frame and is formed with an asymmetrical construction.
- Stirring instruments allow flexion/raising of legs and stirring/string ball joints locked from foot end.
- Two traction devices on sliding spot with ball joint allowing flexions and rotations of foot with traction shoe with metallic sole.
- Padded orthopaedic knee crutches provides support to the leg of the patient.

Standard Accessories

- Boot with block – 1 pair
- Peritoneal post – 1 pc
- Screw tension device with hand crank device – 1 pair
- Hip section with perineal post holder for top traction system – 1 pc
- Cuff support – 1 pair
- Holder for femur lateral to support foot – 1 pc
- Tibia nailing – 1 pc
- Translucent hand rest with adjustable bar – 1 pc
- Steinman pin holder – 1 pc



Tibia Nailing



Femur Nailing



Hip Prosthesis