

Housing	Powder coated body with Stainless Steel Top
Max. Vacuum	-710mm Hg $\pm$ 10mm Hg, regulable through Vacuum Control Switch
Suction Capacity	30 -35 Liters / Min
Vacuum Gauge	2.5 inch , 0 -760mm Hg
Jars	2 x 1.5 Liters Poly Carbonate jars or Glass jars (Optional)
Tubing	Non- Collapsible medical grade tubing
Motor	Motor ¼ hp
Noise Level	<50dB A $\pm$ 3 Almost Silent
Vacuum pump	Rotary vane Type
Certifications	CE Certified

VacuMatrix VM - 30



📍 #487/67, Glass Street, National Market,
Peeragarhi, New Delhi - 110087, INDIA

📞 +91 11 40392484

☎ +91 9812266676

✉ marketinguk@meditech-india.com

🌐 www.mediech-india.com