

Housing	Powder coated body (Trolley Type)
Max. Vacuum	-700mm Hg $\pm$ 10mm Hg, regulable through Vacuum Control Switch
Motor	Oil free double piston Pump
Suction Capacity	35 Liters / Min
Vacuum Gauge	2.5 inch, 0-760 mm Hg
Jars	2 x 1.5 Liters Polycarbonate Jars
Tubing	Non-Collapsible medical grade tubing
Noise Level	<30-35dB A $\pm$ 3 Almost Silent
Overflow Safety	Mechanical
Certifications	CE Certified

VacuMatrix VM - 35



📍 #487/67, Glass Street, National Market,
Peeragarhi, New Delhi - 110087, INDIA

📞 +91 11 40392484

☎ +91 9812266676

✉ marketinguk@meditech-india.com

🌐 www.mediech-india.com