



C-ARM Compatible Spinal Bridge


Salient Features

- Radiolucent and allows for unrestricted C-arm access
- Patient supports are adjustable to fit patients of various sizes
- Adjustable width for adult & child.



Features

- Meditech Spine Frame provides a stable and convenient platform for spine surgeries such as laminectomies, decompressions, and disc surgery
- Patient is positioned in prone position for access to the dorsal surface. Using the included Crank Handle, it can be adjusted intraoperatively to the desired lordosis to get optimal surgical site exposure
- Compatible with any OT Table via loop straps that attach to the table side rail
- Spine Frame: L 26" x W 20" x H 6.5"



MEDITECH (INDIA)[®]
Creating Your Lifeline

487/67, Glass Street, National Market, Peeragarhi, New Delhi-110087, India.
Tel. : +91-11-40392484
Mob. : +91- 9812266676, 9812035477
Email : marketinguk@meditech-india.com
Web : www.meditech-india.com
Skype: Exportmeditech

Design & specifications are subject to change without any further notice. Publishing of any material in full or half of this catalogue is a punishable offence, as the products shown in the catalogue are the sole property of Meditech (India).